



# Orthopaedic Collection Socks

## 20-30 mmHg



### DIAMOND COPPER

EC810C-(BKCO)

20-30 mmHg

COPPER INFUSED FIBERS



### PRIMARY COPPER

EC820C-(COP)

20-30 mmHg

COPPER INFUSED FIBERS



### CITY COPPER

EC860C-(BKCO)

20-30 mmHg

COPPER INFUSED FIBERS




### PRIMARY COPPER

TWIST EC822B-(COP)

20-30 mmHg

COPPER INFUSED FIBERS



### ANKLE COPPER

CG906A-(BKCO)

20-30 mmHg

COPPER INFUSED FIBERS



### TWIST N'GO

EC850C-(BKCO)

20-30 mmHg

COPPER INFUSED FIBERS



### INFINITY

EC830C-(BKW)

20-30 mmHg



### TRENDY

EC840C-(NAB)

20-30 mmHg



### BOREAL

EC842C-(BKG)

20-30 mmHg



### URBAIN

EC844C-(NAP)

20-30 mmHg




#### Measurement Chart for Ankle Socks

		2	3	4	5
<b>N MEN</b>	US [Shoe size]	6-7	7-9	9,5-11,5	12+
<b>P WOMEN</b>	US [Shoe size]	6-8	8-10	10-12	13+

#### Measurement Chart for Compression Socks

<b>K</b> ANKLE MEASUREMENTS	7 in-8,5 in (18-21,5cm)		8,5 in-10 in (21,5-25,5cm)		10 in-11,5 in (25,5-29cm)		11,5 in-14 in (29-36cm)	
<b>H</b> CALF MEASUREMENTS	11 in-15,5 in (28-39cm)		12,5 in-17,5 in (31,1-45cm)		14 in-20 in (35,5-51cm)		16 in-24 in (41-61cm)	
<b>K</b> CALF LENGTH	Short	Long	Short	Long	Short	Long	Short	Long
	< 16 in (41cm)	> 16 in (41cm)	< 16 in (41cm)	> 16 in (41cm)	< 16 in (41cm)	> 16 in (41cm)	< 16 in (41cm)	> 16 in (41cm)
SIZES	2S	2L	3S	3L	4S	4L	5S	5L

#### How to measure





# Collection de Bas Orthopédiques

## 20-30 mmHg



### DIAMOND COPPER

EC810C-(BKCO)

20-30 mmHg

FIBRES INFUSÉES DE CUIVRE



### PRIMARY COPPER

EC820C-(COP)

20-30 mmHg

FIBRES INFUSÉES DE CUIVRE



### CITY COPPER

EC860C-(BKCO)

20-30 mmHg

FIBRES INFUSÉES DE CUIVRE



### PRIMARY COPPER

TWIST EC822B-(COP)

20-30 mmHg

FIBRES INFUSÉES DE CUIVRE



### ANKLE COPPER

CG906A-(BKCO)

20-30 mmHg

FIBRES INFUSÉES DE CUIVRE



### TWIST N'GO

EC850C-(BKCO)

20-30 mmHg

FIBRES INFUSÉES DE CUIVRE



### INFINITY

EC830C-(BKW)

20-30 mmHg



### TRENDY

EC840C-(NAB)

20-30 mmHg



### BOREAL

EC842C-(BKG)

20-30 mmHg



### URBAIN

EC844C-(NAP)

20-30 mmHg



#### Charte de mesure pour Bas Cheville

	2	3	4	5
<b>N HOMME</b> US [Grandeur soulier]	6-7	7-9	9,5-11,5	12+
<b>P FEMME</b> US [Grandeur soulier]	6-8	8-10	10-12	13+

#### Charte de mesure pour les Bas Orthopédiques

<b>O</b> TOUR DE CHEVILLE	7 po-8,5 po (18-21,5cm)		8,5 po-10 po (21,5-25,5cm)		10 po-11,5 po (25,5-29cm)		11,5 po-14 po (29-36cm)	
<b>H</b> TOUR DE MOLLET	11 po-15,5 po (28-39cm)		12,5 po-17,5 po (31,1-45cm)		14 po-20 po (35,5-51cm)		16 po-24 po (41-61cm)	
<b>K</b> LONGUEUR DU MOLLET	Short	Long	Short	Long	Short	Long	Short	Long
	< 16 po (41cm)	> 16 po (41cm)	< 16 po (41cm)	> 16 po (41cm)	< 16 po (41cm)	> 16 po (41cm)	< 16 po (41cm)	> 16 po (41cm)
GRANDEURS	2S	2L	3S	3L	4S	4L	5S	5L

#### Comment Mesurer

